


LITMUS INDICATOR -1




MATERIAL REQUIRED



TAKE SOME FEW GLASS OF LEMON IN WATER GLASS



PUT BLUE LITMUS PAPER ON THE WATCHGLASS THE COLOUR OF LITMUS PAPER CHANGES FROM BLUE TO RED



PUT BLUE LITMUS PAPER ON THE WATCH GLASS, NO CHANGE IN COLOUR OF LITMUS PAPER



TAKE SOME WASHING POWDER DISSOLVED IN WATER IN ANOTHER WATCH GLASS

ARSHA SARAF'S LAB



**PUT RED LITMUS ON THE WATCH GLASS,
COLOUR OF LITMUS PAPER CHANGES
FROM RED TO BLUE**



**PUT RED LITMUS PAPER ON THE WATCH
GLASS, NO CHANGE IN COLOUR OF
LITMUS PAPER**

Litmus paper acts as universal indicator to test acid and bases. In acid, blue litmus turns red. In bases, red litmus turns blue.

ARSHA SARAF'S LAB